

Visceral Leishmaniasis Antibody Detection by Rapid Test

FOR THE QUALITATIVE DETECTION OF ANTI-LEISHMANIA ANTIBODIES

Test Highlights

IgG-specific assay that replaces a total antibodies assay.

Disease Overview

- Visceral leishmaniasis (VL) is also known as kala-azer (Hindi for “black fever”).
- VL is a specific form of leishmaniasis caused by *L. donovani*, *L. chagasi*, and *L. infantum*. Antigenic similarities between these species make IgG testing for *Leishmania donovani* highly sensitive and specific for the disease.
- Parasite is transmitted to humans by sandflies.
- Clinical symptoms of VL are chronic fever, hepatosplenomegaly, anemia, and growth retardation in children; may be fatal if untreated.

Immunology

- Following the onset of symptoms, IgG antibodies are typically measurable in over 90 percent of visceral leishmaniasis cases.
- Testing for IgM antibodies has no demonstrated clinical utility.

Epidemiology

- Visceral leishmaniasis is endemic in the tropics, subtropics, and southern Europe, with the majority of cases occurring in India, Brazil, and Sudan.
- It is estimated that 500,000 people become symptomatic for VL each year, with as many as 350 million people at risk of infection.
- Up to 10 percent of cases remain asymptomatic.
- Military personnel returning from Iraq are at risk for developing VL.
- Infection may be transmitted through blood products.

Indications for Ordering

- This test should be ordered if exposure to sandflies is suspected and visceral leishmaniasis symptoms exist.

- Individuals who have traveled to areas endemic for *Leishmania spp.* may be considered for testing even if asymptomatic, since symptoms may present months following initial infection or not at all.

Interpretation

A positive result of one unit or greater indicates presence of serum IgG antibodies against *Leishmania donovani*.

Limitations

- This test has been validated for serum and plasma samples only. No other sample types may be used.
- This test should not be used in the diagnosis of cutaneous or mucosal leishmaniasis, as an antibody response is not observed in either form of the disease.
- Serological testing has limited utility in acute diagnosis.

Methodology

Patient sample is added to the base of a nitrocellulose strip. The presence of a solid line after incubation indicates a positive result.

References

1. Mandell G, Bennett J, Dolin R, eds. Principles and Practice of Infectious Diseases. 5th ed. 2000; New York: Churchill Livingstone.
2. Davidson R. N. Visceral leishmaniasis in clinical practice. *J Infect* 1999;39:112–6.
3. Herwaldt B. L. Leishmaniasis. *Lancet* 1999;354:1191–9.
4. Package Insert for *Kalazar Detect*™ Rapid Test. InBios International, Inc. 2007.

Test Information

0051726

Leishmania Antibody, IgG (Visceral Leishmaniasis)

For specific collection, transport, and testing information, refer to the ARUP Web site at www.aruplab.com.

For information on test selection, ordering, and interpretation, refer to ARUP Consult® at www.arupconsult.com.