

## Prospera Transplant Assessment (Qualified Clients Only)

### Collection

#### Special Instructions

Required with specimen submission:

- [Patient History Form](#)
- **Patient insurance information**

#### Collection Information

- [Collection kit image](#)
- [Sample Collection Instruction](#)

#### Collect

Whole blood must be collected in 2 Cell-Free DNA BCT tubes included the Prospera Transplant Assessment Collection kit (ARUP Supply #55552) available online through eSupply or by contacting ARUP Client Services at (800) 522-2787.

#### Specimen Preparation

Transport 20 mL blood in Cell-Free DNA BCT Tube (ARUP Supply #55552) available online through eSupply or contacting ARUP Client Services at (800) 522-2787.

#### Storage/Transport Temperature

Room temperature (8 - 37°C)

#### Stability (from collection to initiation)

Room temperature, 5 days

Refrigerated, Unacceptable

Frozen, Unacceptable

#### Unacceptable Conditions

Refrigerated, Frozen

#### Remarks

Do not expose to extreme temperatures

## Prospera Transplant Assessment (Qualified Clients Only)

### Ordering

#### Ordering Recommendations

Assessment of active rejection and injury in kidney transplant recipients that are more than 14 days posttransplant. For outpatient use.

For additional information, please see [Natera's Test Information](#)

#### Methodology

Targeted Sequencing with SNP

#### Performed

Varies

#### Reported

5 - 7 days

#### Synonyms

- kidney transplant
- liquid biopsy
- cfDNA
- dd-cfdNA
- cell-free DNA
- cell free DNA
- PROSPERA

#### Notes (or Technical Information)

Test utilizes a single-nucleotide polymorphism (SNP)/ informatics-based approach to detect donor organ DNA in transplant recipient's blood. The dd-cfDNA fraction is determined using a proprietary algorithm that does not require prior analysis of either donor or recipient DNA. When samples do not meet the necessary quality metrics, a test result is not provided, and the clinician is advised to perform a second draw. False positives and false negatives can occur. High dd-cfDNA fraction, associated with increased risk of active rejection (AR), may require diagnostic confirmation of AR by alternative testing methods. Low dd-cfDNA fraction results do not fully exclude the diagnosis of AR nor do they exclude the possibility of other kidney injuries. Test results should always be interpreted by a clinician in the context of clinical data.

[Additional Technical Information](#)

## Prospera Transplant Assessment (Qualified Clients Only)

### Resulting

#### Reference Interval

By report

#### Interpretive Data

Refer to report

#### Lab Section


Send-Out Lab

### Administrative

#### CPT Codes

N/A (billed by Natera)

#### Referral Lab

 (performed at Natera)

#### Approvals

NY approved

#### Forms & Information

- [Patient History Form](#)
- [Collection information](#)
- [Additional technical information](#)

#### Updated

2-JUL-2020