

A nonprofit enterprise of the University of Utah and its Department of Pathology

## **Order Entry New Label Printing Application**

ARUP is dropping support for the Java applet currently used to print labels in Order Entry and switching to a new printing application. If you use the label printing feature, you will need to download the new application by 5 p.m. on Tuesday, November 3, 2015, and re-setup your label printer.

If you utilize your own LIS labels for specimens that ship to ARUP, please disregard these instructions.

If you do not have administrator rights on your machine, please ask your IT group to assist you with the setup.

We highly recommend that you download the software required for the new label printing application prior to the November 3 release.

Navigate to: http://www.neodynamic.com/downloads/wcpp/

Be sure to download the WCPP-\*.msi file as illustrated below.

ssor (WCPP) is the hines without pay	client-side part of both WebC	lientPrint for ASP.NET & WebClien	tPrint for PHP solutions. It can	n be installed at any
miles manout puy		and read System Requirements for ea	ch available platform:	-
	ing only recsi bounded from t	and read system requirements for ea	ch available platform.	
Vindows	WCPP for Linux	WCPP for Mac		
er File			WCPP Version	
op-2.0.15.918-win	exe		2.0.15.918	
op-2.0.15.918.msi			2.0.15.918	
	<b>Vindows</b> ler File pp-2.0.15.918-win pp-2.0.15.918.msi	Windows WCPP for Linux ler File pp-2.0.15.918-win.exe pp-2.0.15.918.msi	Windows WCPP for Linux WCPP for Mac	Windows         WCPP for Linux         WCPP for Mac           ler File         WCPP Version           pp-2.0.15.918-win.exe         2.0.15.918           pp-2.0.15.918.msi         2.0.15.918

Following our release, you can log into Connect Order Entry and click on the link under the Settings tab to get to the download.

Register Patient	Look Up Patient         Packing List         Packing List Reprint         Label Reprint         Settings
Physicians	Printers
Quick List	✓ Enable Label Printing
Packages	
Printers	Additional software is required for label printing. Please <u>click here</u> to download the software.
	Reboot your computer after installing the application then return to this page and set up your printer.

Double click on the file to run the installation.



Click Next. Select to install for everyone. This is critical if you have multiple users using the same computer.

B WebClientPrint Processor 2.0 for Windows	
Select Installation Folder	
The installer will install WebClientPrint Processor 2.0 for Windows to the following	folder.
To install in this folder, click "Next". To install to a different folder, enter it below or	click "Browse".
Eolder: C:\Program Files (x86)\Neodynamic\WebClientPrint Processor 2.0 for	Browse Disk Cost
Install WebClientPrint Processor 2.0 for Windows for yourself, or for anyone who computer:	) uses this
Cancel < Back	Next >

Click Next to confirm installation.

B WebClientPrint Processor 2.0 for Windows	×
Confirm Installation	5
The installer is ready to install WebClientPrint Processor 2.0 for Windows on you	r computer.
Click "Next" to start the installation.	
Cancel < Back	Next >

Installation will begin.

B WebClientPrint Processor 2.0 for Windows	X
Installing WebClientPrint Processor 2.0 for Windows	-
WebClientPrint Processor 2.0 for Windows is being installed.	
Please wait	
Cancel < Ba	ack Next >

Click Close to finish installation.

B WebClientPrint Processor 2.0 for Windows	X
Installation Complete	5
WebClientPrint Processor 2.0 for Windows has been successfully installed. Click "Close" to exit.	
Cancel < Back	Close

Close all applications and programs and reboot your computer.

Log into Connect Order Entry and select the Settings tab.

Register Patient	Look Up Patient         Packing List         Packing List Reprint         Label Reprint         Settings
Physicians	Printers
Quick List	Inable Label Printing
Packages	Select Printer: Zebra LP2824
Printers	Copies: 1

Use the dropdown menu to select the Zebra LP2824 printer connected to your computer.

You may see the following screen instead of the dropdown menu on some browsers. If you see this screen, simply click on the link to load your printers.

Register Patient	Look Up Patient         Packing List         Packing List Reprint         Label Reprint         Settings
Physicians	Printers
Quick List	C Enable Label Printing
Packages	Could not load the printers.
Printers	Please make sure that printers are installed on your machine. Please <mark>click here</mark> to try loading the installed printers.

Click Save to save your printer selection.

Your computer is now setup for label printing.

## Additional steps for IE 8 only:

Click on the Tools menu and open the Internet Options window as illustrated below.

Select Internet as your zone and press the Custom level button to open the Security Settings – Internet Zone window.



Scroll down to the Miscellaneous section.

Scroll down to the sub-heading "Display mixed content" and select Enable (the default is Prompt).

Press OK, then restart the browser.

Security Settings - Internet Zone	×
Settings	
Prompt	<u> </u>
Allow websites to open windows without address or status b	
Disable	
Enable	
Display mixed content	
Disable	
Enable	
Prompt	
Don't prompt for client certificate selection when no certificat	
Oisable	
Enable	
Drag and drop or copy and paste files	
Disable	
<ul> <li>Enable</li> </ul>	
Prompt	-
Todude local directory path when uploading files to a server	·
*Takes effect after you restart Internet Explorer	
Reset custom settings	
Reset to: Medium high (default)	
reduin nigh (derdale)	
OK	el 🚽