



Quarterly HOT LINE: Effective **August 15, 2016**

0092420	Drug Screen 9 Panel, Serum or Plasma - Immunoassay Screen with Reflex to Mass Spectrometry Confirmation/Quantitation	DRUG SCRSP
----------------	---	-------------------

Reference Interval:

Drugs Covered and Cutoff Concentrations

Drugs/Drug Classes	Screen
Amphetamines	30 ng/mL
Methamphetamine	30 ng/mL
Barbiturates	75 ng/mL
Benzodiazepines	75 ng/mL
Buprenorphine	Effective Aug 15, 2016 1 ng/mL
Cannabinoids	30 ng/mL
Cocaine	30 ng/mL
Methadone	40 ng/mL
Opiates	30 ng/mL
Oxycodone	30 ng/mL
Phencyclidine	15 ng/mL

CPT Code(s): 80301; if positive, add appropriate CPT code(s): 80324; 80359; 80345; 80347; 80349; 80353; 80358; 80361; 80365; **80348**; 83992 (Alt code: G0479; if positive, add appropriate CPT code(s): G0480)

HOT LINE NOTE: There is a component change associated with this test.

Remove component 0092429, Propoxyphene, S/P, Screen

Add component 2013568, Buprenorphine, S/P, Screen