Example Report

ARUP Laboratories
500 Chipeta Way – Salt Lake City, UT 84108

(800)522-2787 - www.aruplab.com Julio C. Delgado, M.D. M.S., Director of Laboratories Patient Age/Gender: 33 years Female Printed: 23-Jun-17 11:00:29

 Procedure
 Result Candida Antibody IgG
 EV
 Inits EV
 Ref Interval EV
 Accession Interval EV
 Accession Interval Inter

14-Jun-17 15:15:00 Candida Antibody IgG: INTERPRETIVE INFORMATION: Candida Ab, IgG

0.89 EV or less: Negative - No significant level

of detectable Candida albicans

antibody.

0.90-0.99 EV: Equivocal - Questionable presence

of antibodies. Repeat testing in

10-14 days may be helpful.

1.00 EV or greater: Positive - Antibody to Candida

albicans detected, which may indicate a current or past

infection.

The best evidence for current infection is a significant change on two appropriately timed specimens where both tests are done in the same laboratory at the same time. However, low levels of IgM antibodies may occasionally persist for more than 12 months post-infection.

Test developed and characteristics determined by ARUP Laboratories. See Compliance Statement D: aruplab.com/CS

Chart ID: 12327697 Page 1 of 1