Example Report

ARUP Laboratories 500 Chipeta Way - Salt Lake City, UT 84108 (800)522-2787 - www.aruplab.com Julio C. Delgado, M.D. M.S., Director of Laboratories

Patient Age/Gender: 19 years Male Printed: 20-Mar-19 10:48:45

Procedure Result Ref Interval Accession Units $\frac{\text{Accession}}{19\text{-}079\text{-}900092} \, \frac{\text{Collected}}{20\text{-}\text{Mar}\text{-}19} \, \frac{\text{Received}}{20\text{-}\text{Mar}\text{-}19} \, \frac{\text{Verified}}{20\text{-}\text{Mar}\text{-}19}$ HIV-1 Qualitative by NAAT Detected Not

Detected]

10:42:00 10:43:00 10:43:52

20-Mar-19 10:42:00 HIV-1 Qualitative by NAAT: INTERPRETIVE DATA: HIV-1 by Qualitative NAAT, Whole Blood

This test detects human immunodeficiency virus type 1 (HIV-1) RNA from Group M, N and O subtypes.

A result of "Not Detected" does not rule out HIV-1 RNA concentrations below the limit of detection of the assay or the presence of inhibitors in the patient specimen. This assay may not detect HIV infection in infants during the first months of life.

The diagnosis of HIV-1 infection should be made based on clinical presentation and results from additional diagnostic tests. Diagnosis should not be made based solely on a single HIV-1 test.

Improper specimen handling can cause false negatives or contamination.

This assay should not be used for blood donor screening, associated re-entry protocols, or for screening Human Cell, Tissues and Cellular Tissue-Based Products (HCT/P).

Test developed and characteristics determined by ARUP Laboratories. See Compliance Statement B: aruplab.com/CS

* Abnormal, # = Corrected, C = Critical, f = Footnote, H = High, L = Low, t = Interpretive Text, @ = Reference Lab

Chart ID: 13330148 Page 1 of 1