

Effective as of **10/05/2026**

Additional ordering and billing information

[Information when ordering laboratory tests that are billed to Medicare/Medicaid](#)

[Information regarding Current Procedural Terminology \(CPT\)](#)

Test Number	Mnemonic	Test Name	New Test	Test Name Change	Specimen Requirements	Methodology	Note	Interpretive Data	Reference Interval	Component Charting Name	Component Change	Reflex Pattern	Result Type	Ask at Order Prompt	Numeric Map	Unit of Measure	CPT Code	Pricing Change	Inactivation w/ Replacement	Inactivation w/o Replacement
2005593	NATAL ABS	Natalizumab Antibodies																		x
2006593	PORPH TOT	Porphyrins, Total with Reflex to Porphyrins Fractionation, Plasma			x															
3002337	BETA PG U	2,3 Dinor-11Beta-Prostaglandin F2 Alpha, 24-Hour Urine			x															
3004160	BETAPG RAN	2,3 Dinor-11Beta-Prostaglandin F2 Alpha, Random Urine			x															

TEST CHANGE

Porphyrins, Total with Reflex to Porphyrins Fractionation, Plasma

2006593, PORPH TOT

Specimen Requirements:

Patient Preparation: Patient should abstain from alcohol for 24 hours prior to collection.

Collect: Green (sodium or lithium heparin). Also acceptable: lavender (EDTA).

Specimen Preparation: Protect from light. Transfer 3 mL plasma to an ARUP amber transport tube ~~and freeze immediately~~. (Min: 1 mL)
Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

Transport Temperature: ~~CRITICAL FROZEN~~. Frozen. Also acceptable: Refrigerated

Unacceptable Conditions: Separator tubes. Specimens not protected from light.

Remarks: Include a current list of patient medications.

Stability: Ambient: Unacceptable; Refrigerated: ~~Unacceptable~~ 2 weeks; Frozen: 2 weeks

Methodology: ~~Qualitative Extraction /~~ Quantitative Scanning Spectrofluorometry / Quantitative High Performance Liquid Chromatography (HPLC)

Note: If total porphyrins are greater than 1.0 mcg/dL, then fractionation will be performed. Additional charges apply.

CPT Codes: 84311; if reflexed, add 82542

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Reference Interval:

~~By Report~~ Refer to report

TEST CHANGE

2,3 Dinor-11Beta-Prostaglandin F2 Alpha, 24-Hour Urine

3002337, BETA PG U

Specimen Requirements:

Patient Preparation: Patients taking aspirin or nonsteroidal anti-inflammatory drugs (NSAIDs) may have decreased concentrations of prostaglandin F2 alpha. If ~~possible, discontinue~~medically feasible, for 2 weeks ~~or 72 hours, respectively~~, prior to specimen collection. patient should not take aspirin and for 72 hours prior to specimen collection, patient should not take NSAIDs.

Collect: 24-hour urine. Refrigerate during collection.

Specimen Preparation: From a well-mixed 24-hour collection transfer 5 mL urine to ARUP standard transport tubes. (Min: ~~2-73~~ mL)
Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

Transport Temperature: Refrigerated. Also acceptable: Frozen.

Unacceptable Conditions:

Remarks: Collection duration and urine volume must be provided for testing. Total collection volume must be greater than 300 mL.

Stability: Ambient: ~~8 hours~~3 days; Refrigerated: ~~2 weeks~~28 days; Frozen: ~~1 month~~28 days

Methodology: Quantitative Colorimetry / ~~High Performance~~Quantitative Liquid Chromatography-Tandem Mass Spectrometry

Note: Elevated levels of 2,3-dinor-11beta-prostaglandin F2 alpha (2,3 BPG) in urine are not specific for systemic mast cell disease and may be found in patients with angioedema, diffuse urticaria, or myeloproliferative diseases in the absence of diffuse mast cell proliferation.

CPT Codes: 84150

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Reference Interval:

~~By Report~~Refer to report

TEST CHANGE

2,3 Dinor-11Beta-Prostaglandin F2 Alpha, Random Urine

3004160, BETAPG RAN

Specimen Requirements:

Patient Preparation: Patients taking aspirin or nonsteroidal anti-inflammatory drugs (NSAIDs) may have decreased concentrations of prostaglandin F2 alpha. If ~~possible, discontinued~~medically feasible, for 2 weeks ~~or 72 hours, respectively~~, prior to specimen collection. patient should not take aspirin and for 72 hours prior to specimen collection, patient should not take NSAIDs.

Collect: Urine.

Specimen Preparation: Transfer 5 mL urine to ARUP standard transport tubes. (Min: 3 mL)
Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

Transport Temperature: Refrigerated. Also acceptable: Frozen.

Unacceptable Conditions:

Remarks:

Stability: Ambient: ~~8 hours~~3 days; Refrigerated: ~~2 weeks~~28 days; Frozen: ~~1 month~~28 days

Methodology: Quantitative Liquid Chromatography-Tandem Mass Spectrometry / Colorimetry

Note: Elevated levels of 2,3-dinor-11beta-prostaglandin F2 alpha (2,3 BPG) in urine are not specific for systemic mast cell disease and may be found in patients with angioedema, diffuse urticaria, or myeloproliferative diseases in the absence of diffuse mast cell proliferation.

CPT Codes: 84150; 82570

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Reference Interval:

~~By Report~~Refer to report

Inactivations

The following will be discontinued from ARUP's test menu on **October 5, 2026**

Replacement test options are indicated when applicable.

Test Number	Test Name	Refer to Replacement Test
2005593	Natalizumab Antibodies	