

IMMEDIATE CHANGE HOTLINE: Effective May 1, 2017



Delete

2007780

Allergen, Food, Cumin Seed IgE

CUMIN

THROMBO S

HOTLINE NOTE: Delete this test.

2013335 Thrombopoietin (TPO), Serum

Performed:VariesReported:7-40 days

Specimen Required: Collect: Plain Red. Also acceptable: Serum Separator Tube (SST).

<u>Specimen Preparation:</u> Transfer 1 mL serum to an ARUP Standard Transport Tube. (Min: 0.5 mL) <u>Storage/Transport Temperature:</u> **CRITICAL FROZEN. Separate specimens must be submitted when multiple tests are ordered.** <u>Unacceptable Conditions:</u> Thawed specimens. <u>Stability (collection to initiation of testing):</u> Ambient: Unacceptable; Refrigerated: Unacceptable; Frozen: 1 month