

Laboratory Diagnosis of Hemoglobinopathies and Thalassemia

Archana M Agarwal, MD

Medical Director, Hematopathology and RBC
Laboratory
ARUP Laboratories
Assistant Professor of Pathology
University of Utah Department of Pathology

ARUP LABORATORIES | NATIONAL REFERENCE LABORATORY

UNIVERSITY OF UTAH
HARVEY C. HARRIS
Department of Pathology

Learning Objectives

- Understand the pathophysiology of hemoglobinopathies
- Recognize the most important expected test results in hemoglobinopathies and thalassemias
- Understand different testing methodologies
- To be able to direct ordering physician to appropriate tests for these disorders

ARUP LABORATORIES | NATIONAL REFERENCE LABORATORY

UNIVERSITY OF UTAH
HARVEY C. HARRIS
Department of Pathology


Hemoglobin (Heme+Globin)

- Hemoglobin is a tetramer composed of 4 globin molecules; 2 alpha globins and 2 beta globins or beta like globins
- The alpha globin chain is composed of 141 amino acids and the beta globin chain is composed of 146 amino acids
- Each globin chain also contains one heme molecule

ARUP LABORATORIES | NATIONAL REFERENCE LABORATORY

UNIVERSITY OF UTAH
HARVEY C. HARRIS
Department of Pathology

Ribbon Diagram of Hemoglobin




ARUP

NATIONAL REFERENCE LABORATORY


UNIVERSITY OF UTAH
HARVEY C. HARRIS SCHOOL OF MEDICINE
Department of Pathology

Genetics of Globin Genes

Chromosome 16



Chromosome 11




ARUP

NATIONAL REFERENCE LABORATORY

UNIVERSITY OF UTAH
HARVEY C. HARRIS SCHOOL OF MEDICINE
Department of Pathology

Hemoglobin-Development Switching



ARUP

NATIONAL REFERENCE LABORATORY

UNIVERSITY OF UTAH
HARVEY C. HARRIS SCHOOL OF MEDICINE
Department of Pathology

Normal Adult Human Hemoglobin Composition

Hemoglobin	Structure	% of Normal Adult Hb
Hb A	$\alpha_2\beta_2$	>96%
Hb A2	$\alpha_2\delta_2$	~2.5%
Hb F	$\alpha_2\gamma_2$	<1%

Hemoglobinopathy (structural)

- Due to mutations in either alpha or beta globin
- **Structural** – substitution, addition or deletion of one or more AAs in the globin chain
 - i.e HbS, HbC, HbE, HbD, HbO, etc...
- Over 1000 identified
 - Majority are benign & discovered incidentally
 - Pathogenic mutations can cause
 - Change in physical properties (sickling, crystalizes)
 - Globin instability (Heinz body formation, lower expression)
 - Altered oxygen affinity

Thalassemia (quantitative)

- A quantitative decrease in the production of alpha or beta globin chain
 - Large deletions, point mutations, small insertion/deletion that leads to decreased transcription or an unstable transcript
- Beta thalassemia results from mutations in beta gene(s)
 - Pathogenesis a result of the **free alpha subunits**
 - Two classes: **β^0 and β^+**
- Alpha thalassemia results from large deletions in the alpha gene(s)
 - Pathogenesis a result of the **free beta subunits**

Demographics: Thalassemias

- Found most frequently in the Mediterranean, Africa, Western and Southeast Asia, India and Burma
- Distribution parallels that of *Plasmodium falciparum*



ARUP

NATIONAL REFERENCE LABORATORY

UNIVERSITY OF UTAH
HARVEY C. HARRIS
Department of Pathology

Classification & Terminology: Alpha Thalassemia

- Normal $\alpha\alpha/\alpha\alpha$
- Silent carrier $-\alpha/\alpha\alpha$
- Minor $-\alpha/-\alpha$ trait $--/\alpha\alpha$
- Hb H disease $--/-\alpha$
- Barts hydrops fetalis $--/--$

ARUP

NATIONAL REFERENCE LABORATORY

UNIVERSITY OF UTAH
HARVEY C. HARRIS
Department of Pathology

Clinical Presentations of Alpha Thalassemia

- A **single** deletion (α -thalassemia minor)
 - silent carrier state
 - RBC morphology and hemoglobin concentrations are usually normal
- **Two** gene deletion (α -thalassemia minor)
 - Mild microcytic anemia
- **Three** gene deletion (**hemoglobin H disease**)
 - Precipitated β chains—Hb H
 - Patients have moderate anemia, marked microcytosis, splenomegaly, and bone marrow erythroid hyperplasia
- **Four** gene deletion (Hydrops fetalis)
 - Not compatible with life (barring very early intervention)
 - Hemoglobin is primarily comprised of $\gamma 4$ (Bart's), which has a very high affinity for O₂ and is a poor oxygen transporter

ARUP

NATIONAL REFERENCE LABORATORY

UNIVERSITY OF UTAH
HARVEY C. HARRIS
Department of Pathology

Classification & Terminology:

Beta Thalassemia

- Normal β/β
- Minor / trait β/β^0
 β/β^+
- Intermedia β^0/β^+
- Major β^0/β^0
 β^+/ β^+


Clinical Significance of β Thalassemia

- Heterozygous asymptomatic
- Homozygous β^0 is a severe disorder associated with transfusion dependent hemolytic anemia
- Homozygous β^+ is a heterogenous disorder
 - severity depending on mutation and % of HbA
 - Increased HbA = decreased severity

Sickle Cell Anemia

- Single nucleotide base change codes for valine instead of glutamic acid at the 6th position from the N-terminus of the β -globin chain
- Affects the shape and deformability of the red blood cell
- Leads to veno-occlusive disease and hemolysis

Peripheral Smear: Sickle Cell Anemia



ARUP LABORATORIES
NATIONAL REFERENCE LABORATORY

UNIVERSITY OF UTAH
HARVEY M. GARDNER
Department of Pathology

Hb E


- 2nd most prevalent hemoglobin variant
 - 30,000,000 world wide
 - 80% in Southeast Asia
- Hb E trait: microcytosis (mean MCV=65fl). No anemia
- Hb E disease: MCV =55-65fl with minimal anemia
- *On HPLC has similar migration pattern as Hb A2

ARUP LABORATORIES
NATIONAL REFERENCE LABORATORY

UNIVERSITY OF UTAH
HARVEY M. GARDNER
Department of Pathology

Hb C

- Mutation in β -globin gene $\beta(6\text{glu} \rightarrow \text{lys})$
- Seen predominantly in blacks: Gene prevalence in US black population is 2 to 3%
- May confer malaria resistance
- Often asymptomatic, mild anemia, splenomegaly
- Blood smear shows many target cells, rare intracellular crystals
- Hb S/C disease causes moderate to severe anemia and hemolysis



ARUP LABORATORIES
NATIONAL REFERENCE LABORATORY

UNIVERSITY OF UTAH
HARVEY M. GARDNER
Department of Pathology

Diagnosis

- **Indications for Testing**

- Hemolytic anemia; family history of hemoglobinopathy

- **Laboratory Testing**

- Initial testing – CBC with peripheral smear
- Polychromasia, spherocytes, schistocytes, sickle cells, Heinz bodies, basophilic stippling; however, the lack of any of these cells does not rule out hemolytic anemia
- Many hemoglobinopathies can be diagnosed using electrophoretic or high performance liquid chromatography (HPLC) techniques, but some may be missed
- Genetic testing

Importance of CBC


- **Thalassemias**

- Red cell indices are critical to diagnosis
- Hypochromic microcytic anemia
 - MCV (mean corpuscular volume or size of the cell) is key
 - RDW (red cell distribution width) changes are variable
 - Increased RBC count → one distinguishing factor between thalassemias and other microcytic anemias

Distinguishing Features Between Iron Deficiency and Thalassemia

- | | |
|---|-------------------------------|
| • The RBC count in thalassemia is either normal or on higher side of normal | • Low RBC count |
| • MCV usually less than 70 in | • MCV usually more than 70 |
| • The RDW is usually in the normal range | • RDW is usually more than 17 |

Diagnosis of Thalassemias



ARUP

NATIONAL REFERENCE LABORATORY

UNIVERSITY OF UTAH
HARVEY M. HARRIS
Department of Pathology

High-Pressure Liquid Chromatography



- Cation Exchange
- Analytical cartridge contains negatively charged silica
- Buffers contain Na⁺ and K⁺ ions
- Hemolysates contain positively charged hemoglobin
- Hemoglobin binds to negatively charged silica at injection
- Na⁺ and K⁺ concentration increased and separates hemoglobin fragments from silica

ARUP

NATIONAL REFERENCE LABORATORY

UNIVERSITY OF UTAH
HARVEY M. HARRIS
Department of Pathology

Normal Patient Chromatograms



ARUP

NATIONAL REFERENCE LABORATORY

UNIVERSITY OF UTAH
HARVEY M. HARRIS
Department of Pathology

Summary of HPLC



Advantages

- Fast
- Small amounts of sample
- Accurate quantitation of A2

Disadvantages


- Hemoglobin E cannot be separated from A2
- Hemoglobin H and Barts elute too quickly from column

Capillary Electrophoresis




<http://www.sebia-usa.com>

Phoresis Reports




<http://www.sebia-usa.com>



<http://www.sebia-usa.com>

Alkaline and Acid Gel Electrophoresis

- Electrophoresis (pH 8.4 (alkaline) and pH 6.2 (acid) on agarose gels)
- Slow, labor-intensive, and inaccurate in the quantification of low-concentration Hb variants (e.g., Hb A₂) or in the detection of fast Hb variants (Hb H, Hb Barts)
- The precision and accuracy of Hb A₂ measurements using densitometric scanning of electrophoretic gels is poor, especially when compared with HPLC techniques



Isoelectric Focusing

- IEF is an electrophoretic technique with excellent resolution
- IEF is an equilibrium process in which Hb migrates in a pH gradient to a position of 0 net charge
- The Hb migration order of IEF is the same as that of **alkaline electrophoresis** with better resolution

Molecular Analysis

- Alpha thalassemia
 - Multiplex ligation dependent probe amplification (MLPA) and multiplex PCR
 - Alpha globin sequencing
- Beta thalassemia
 - Beta globin sequencing
 - The test examines the complete beta globin coding sequence, the splice sites and other intronic regions known to harbor mutations, the proximal promoter region, and the 5' and 3'UTR regions.
 - Clinical sensitivity is up to 97% based on the ethnicity
 - Beta globin del/dup testing by MLPA

α -Thalassemia Diagnosis

- Hb gel/HPLC migration patterns
 - Not helpful for α -Thalassemia, unless $\beta 4$ (Hb H) and $\gamma 4$ (Hb Barts) are present
- Genetic analysis
 - MLPA: will identify all deletions and duplications
 - Multiplex PCR for 7 common deletions-only 7 common deletion
 - Alpha globin sequencing
 - PCR amplification followed by bidirectional sequencing of the complete protein coding sequence with exon/intron boundaries, proximal promoter region, 5' and 3' untranslated regions, and polyadenylation signal
 - Only useful in 5-10% of cases where alpha thal is due to point mutation

β -Thalassemia Diagnosis


- **HPLC**: Elevated HB A2 diagnostic
- **Molecular analysis**: Complete beta globin coding sequence, the splice sites and other intronic regions known to harbor mutations, the proximal promoter region, and the 5' and 3'UTR regions
- Clinical sensitivity is up to 97% based on the ethnicity
- Beta globin del/dup in some cases (about 5%) where beta thalassemia is due to large deletions

Sickle Cell Disease Diagnosis

- Sickledex test (Screening test)
 - Deoxygenated Hb-S is insoluble in a concentrated phosphate buffer solution and forms a turbid suspension
 - Normal Hemoglobin A and other hemoglobins remain in solution
 - It does not differentiate between Sickle Cell Disease (S/S) and Sickle Cell Trait (A/S)




Sickle Cell Disease Diagnosis



Color Atlas of Hemoglobin Disorders: A compendium Based on Proficiency Testing (2003), updated in 2010

Simplified Algorithm



References and Acknowledgement

- Kaushansky K, Lichtman MA, Beutler E, Kipps TJ, Prchal J, Seligsohn U. William's Hematology. Ninth Edition. McGraw Hill Professional. 2015.
- Steinberg MH, Forget BG, Higgs DR, Nagel RL. Disorders of Hemoglobin. Genetics, Pathophysiology, and Clinical Management, 2nd ed. Cambridge University Press, New York, 2009
- Color Atlas of Hemoglobin Disorders: A compendium Based on Proficiency Testing (2003), updated in 2010
- **Acknowledgement:**
 - Josef T. Prchal, M.D, Professor of Medicine, Genetics and Pathology. University of Utah and ARUP Laboratories
 - Dottie Hussie, M.T, ARUP Laboratories



© 2014 ARUP Laboratories

ARUP IS A NONPROFIT ENTERPRISE OF THE UNIVERSITY OF UTAH AND ITS DEPARTMENT OF PATHOLOGY.
