



<u>Procedure</u>	<u>Result</u>	<u>Units</u>	<u>Ref Interval</u>	<u>Accession</u>	<u>Collected</u>	<u>Received</u>	<u>Reported/</u> <u>Verified</u>
Total Protein, Serum	7.9 H	g/dL	[5.9-7.8]	19-158-900102	07-Jun-19 10:06:00	07-Jun-19 10:06:00	07-Jun-19 11:59:36
Albumin	4.00	g/dL	[3.75-5.01]	19-158-900102	07-Jun-19 10:06:00	07-Jun-19 10:06:00	07-Jun-19 11:59:36
Alpha 1 Globulin	0.25	g/dL	[0.19-0.46]	19-158-900102	07-Jun-19 10:06:00	07-Jun-19 10:06:00	07-Jun-19 11:59:36
Alpha 2 Globulin	0.50	g/dL	[0.48-1.05]	19-158-900102	07-Jun-19 10:06:00	07-Jun-19 10:06:00	07-Jun-19 11:59:36
Beta Globulin	0.50	g/dL	[0.48-1.10]	19-158-900102	07-Jun-19 10:06:00	07-Jun-19 10:06:00	07-Jun-19 11:59:36
Gamma	0.63	g/dL	[0.62-1.51]	19-158-900102	07-Jun-19 10:06:00	07-Jun-19 10:06:00	07-Jun-19 11:59:36
Immunofixation Reflex	IFE Done			19-158-900102	07-Jun-19 10:06:00	07-Jun-19 10:06:00	07-Jun-19 11:59:36
SPEP/IFE Interpretation	See Note			19-158-900102	07-Jun-19 10:06:00	07-Jun-19 10:06:00	07-Jun-19 11:59:36
EER Serum Protein Electrophoresis Reflex	See Note f			19-158-900102	07-Jun-19 10:06:00	07-Jun-19 10:06:00	07-Jun-19 11:59:36

07-Jun-19 10:06:00 EER Serum Protein Electrophoresis Reflex:
 Access ARUP Enhanced Report using either link below:

-Direct access:

-Enter Username, Password: <https://c11-erpt.aruplab.com>
 Username:
 Password:

* Abnormal, # = Corrected, C = Critical, f = Footnote, H = High, L = Low, t = Interpretive Text, @ = Reference Lab