

**Specimen Collected: 30-Nov-21 11:30**

Amylase Isoenzymes		Received: 30-Nov-21 11:30	Report/Verified: 30-Nov-21 11:33
Procedure	Result	Units	Reference Interval
Pancreatic Amylase Isoenzyme	54 <sup>H</sup> <sup>i1</sup>	U/L	13-53
Salivary Amylase Isoenzyme	47 <sup>i2</sup>	U/L	9-86
Total Amylase for Isoenzymes	101 <sup>H</sup> <sup>i3</sup>	U/L	28-100

**Test Information**

i1: Pancreatic Amylase Isoenzyme  
 REFERENCE INTERVAL: Pancreatic Amylase Isoenzyme

Access complete set of age- and/or gender-specific reference intervals for this test in the ARUP Laboratory Test Directory (aruplab.com).

i2: Salivary Amylase Isoenzyme  
 REFERENCE INTERVAL: Salivary Amylase Isoenzyme

Access complete set of age- and/or gender-specific reference intervals for this test in the ARUP Laboratory Test Directory (aruplab.com).

i3: Total Amylase for Isoenzymes  
 REFERENCE INTERVAL: Total Amylase for Isoenzymes

Access complete set of age- and/or gender-specific reference intervals for this test in the ARUP Laboratory Test Directory (aruplab.com).

\*=Abnormal, #=Corrected, C=Critical, f=Result Footnote, H=High, i=Test Information, L=Low, t=Interpretive Text, @=Performing lab

**Unless otherwise indicated, testing performed at:**

**ARUP Laboratories**  
 500 Chipeta Way, Salt Lake City, UT 84108  
 Laboratory Director: Tracy I. George, MD

**ARUP Accession:** 21-334-900103  
**Report Request ID:** 15063635  
**Printed:** 07-Dec-21 15:34  
 Page 1 of 1