



| <u>Procedure</u> | <u>Result</u> | <u>Units</u> | <u>Ref Interval</u> | <u>Accession</u> | <u>Collected</u>   | <u>Received</u>    | <u>Reported/</u><br><u>Verified</u> |
|------------------|---------------|--------------|---------------------|------------------|--------------------|--------------------|-------------------------------------|
| Arsenic Blood    | >500.0 H      | ug/L         | [0.0-12.0]          | 18-254-900046    | 11-Sep-18 08:15:00 | 11-Sep-18 08:15:00 | 11-Sep-18 10:01:53                  |

11-Sep-18 08:15:00 Arsenic Blood:  
 INTERPRETIVE INFORMATION: Arsenic, Blood

Elevated results may be due to skin or collection-related contamination, including the use of a noncertified metal-free collection/transport tube. If contamination concerns exist due to elevated levels of blood arsenic, confirmation with a second specimen collected in a certified metal-free tube is recommended.

Potentially toxic ranges for blood arsenic: Greater than or equal to 600 ug/L.

Blood arsenic is for the detection of recent exposure poisoning only. Blood arsenic levels in healthy subjects vary considerably with exposure to arsenic in the diet and the environment. A 24-hour urine arsenic is useful for the detection of chronic exposure.

Test developed and characteristics determined by ARUP Laboratories. See Compliance Statement B: aruplab.com/CS

\* Abnormal, # = Corrected, C = Critical, f = Footnote, H = High, L = Low, t = Interpretive Text, @ = Reference Lab