

Specimen Collected: 16-Sep-21 11:35

Heavy Metals Panel 3, Urine with Reflex | Received: 16-Sep-21 11:35 Report/Verified: 16-Sep-21 11:51

Procedure	Result	Units	Reference Interval
Hours Collected	24 <sup>f1</sup>	hr	
Total Volume	1000	mL	
Creatinine, Urine -per volume	100	mg/dL	
Creatinine, Urine -per 24h	1000	mg/d	700-1600
Lead, Urine -per volume	<5.0 <sup>i1</sup>	ug/L	0.0-5.0
Lead, Urine -per 24h	Not Applicable	ug/d	0.0-8.1
Lead, Urine -ratio to CRT	Not Applicable <sup>f2</sup>	ug/g CRT	0.0-5.0
Mercury, Urine -per volume	<2.5 <sup>i2</sup>	ug/L	0.0-5.0
Mercury, Urine -per 24h	Not Applicable	ug/d	0.0-20.0
Mercury, Urine -ratio to CRT	Not Applicable <sup>f2</sup>	ug/g CRT	0.0-20.0
Arsenic Urine -per volume	50.0 <sup>H i3</sup>	ug/L	0.0-34.9
Arsenic Urine -per 24h	50.0 <sup>H</sup>	ug/d	0.0-49.9
Arsenic, Urine -ratio to CRT	50.0 <sup>H</sup>	ug/g CRT	0.0-29.9

Procedure	Result	Units	Reference Interval
Arsenic, Organized, Urine	50.0	ug/L	
Arsenic, Inorganic	<10.0	ug/L	
Arsenic, Methylated	<10.0 <sup>i4</sup>	ug/L	

**Result Footnote**

f1: Hours Collected

Per 24h calculations are provided to aid interpretation for collections with a duration of 24 hours and an average daily urine volume. For specimens with notable deviations in collection time or volume, ratios of analytes to a corresponding urine creatinine concentration may assist in result interpretation.

f2: Lead, Urine - ratio to CRT, Mercury, Urine - ratio to CRT

Unable to accurately calculate the creatinine normalized result due to a low per volume result.

**Test Information**

i1: Lead, Urine - per volume

INTERPRETIVE INFORMATION: Lead, Urine

Quantification of urine excretion rates before or after chelation therapy has been used as an indicator of lead exposure. Urinary excretion of >125 mg of lead per 24 hours is usually associated with related evidence of lead toxicity.

\*=Abnormal, #=Corrected, C=Critical, f=Result Footnote, H=High, i=Test Information, L=Low, t=Interpretive Text, @=Performing lab

Unless otherwise indicated, testing performed at:

ARUP Laboratories

500 Chipeta Way, Salt Lake City, UT 84108

Laboratory Director: Tracy I. George, MD

ARUP Accession: 21-259-900063

Report Request ID: 15048343

Printed: 16-Sep-21 12:24

Page 1 of 3

**Test Information**

i1: Lead, Urine - per volume

This test was developed and its performance characteristics determined by ARUP Laboratories. It has not been cleared or approved by the US Food and Drug Administration. This test was performed in a CLIA certified laboratory and is intended for clinical purposes.

i2: Mercury, Urine - per volume

INTERPRETIVE INFORMATION: Mercury, Urine

Urinary mercury levels predominantly reflect acute or chronic elemental or inorganic mercury exposure. Urine concentrations in unexposed individuals are typically less than 10 ug/L. 24 hour urine concentrations of 30 to 100 ug/L may be associated with subclinical neuropsychiatric symptoms and tremors. Concentrations greater than 100 ug/L can be associated with overt neuropsychiatric disturbances and tremors. Urine mercury levels may be useful in monitoring chelation therapy.

This test was developed and its performance characteristics determined by ARUP Laboratories. It has not been cleared or approved by the US Food and Drug Administration. This test was performed in a CLIA certified laboratory and is intended for clinical purposes.

i3: Arsenic Urine - per volume

INTERPRETIVE INFORMATION: Arsenic, Urine w/ Reflex to Fractionated

The ACGIH Biological Exposure Index (BEI) for arsenic in urine is 35 µg/L. The ACGIH BEI is based on the sum of inorganic and methylated species. For specimens with elevated total arsenic results, fractionation is automatically performed to determine the proportions of inorganic, methylated and organic species.

This test was developed and its performance characteristics determined by ARUP Laboratories. It has not been cleared or approved by the US Food and Drug Administration. This test was performed in a CLIA certified laboratory and is intended for clinical purposes.

i4: Arsenic, Methylated

INTERPRETIVE INFORMATION: Arsenic, Fractionated Urine

The ACGIH Biological Exposure Index for the sum of inorganic and methylated species of arsenic is 35 ug/L. Inorganic species of arsenic are most toxic. Methylated species arise primarily from metabolism of inorganic species but may also come from dietary sources and are of moderate toxic potential. The organic species of arsenic are considered nontoxic and arise primarily from food. The sum of the inorganic, methylated, and organic species of arsenic may be lower than the total arsenic concentration due to the presence of unidentified organic species of arsenic.

This test was developed and its performance characteristics determined by ARUP Laboratories. It has not been cleared or approved by the US Food and Drug

\*=Abnormal, #=Corrected, C=Critical, f=Result Footnote, H=High, i=Test Information, L=Low, t=Interpretive Text, @=Performing lab

**Unless otherwise indicated, testing performed at:**

**ARUP Laboratories**

500 Chipeta Way, Salt Lake City, UT 84108

Laboratory Director: Tracy I. George, MD

**ARUP Accession:** 21-259-900063

**Report Request ID:** 15048343

**Printed:** 16-Sep-21 12:24

Page 2 of 3

**Test Information**

i4: Arsenic, Methylated

Administration. This test was performed in a CLIA certified laboratory and is intended for clinical purposes.

---

\*=Abnormal, #=Corrected, C=Critical, f=Result Footnote, H=High, i=Test Information, L=Low, t=Interpretive Text, @=Performing lab

---

***Unless otherwise indicated, testing performed at:***

**ARUP Laboratories**

500 Chipeta Way, Salt Lake City, UT 84108

Laboratory Director: Tracy I. George, MD

**ARUP Accession:** 21-259-900063

**Report Request ID:** 15048343

**Printed:** 16-Sep-21 12:24

Page 3 of 3