

CLIENT SUPPLY CHECKLIST

Date: ____/____/____ Client #: _____

Ship to address: _____

Ordered by: _____
(full name)

Attn: _____

Phone #: _____

FAX to 1 (800) 522-2706.

Sections in RED are REQUIRED. Missing information will
delay order processing.

email: _____

Please allow five to seven business days for delivery.

QUANTITY

BAGS

_____	Bag, Sample Amber Ziplock PK/100 (4x6)	10528
_____	Bag, Sample Frozen PK/50 (Blue)	22274
_____	Bag, Sample Refrigerated PK/50 (Grey)	22273
_____	Bag, Sample Room Temperature PK/50 (Green)	22272
_____	Bag, Sample Biohazard re-sealable PK/100 6x6	42315

CYTOLOGY SUPPLIES

_____	Broom, Cervical PK/100	22218
_____	Brush, Cervical	11440
_____	Spatula, Cervical PK/100	22219
_____	Kit, SurePath™ 1 kit = 25 vials/Rovers brushes #	22216
_____	Kit, SurePath™ 1 kit = 25 vials/brushes/spatulas #	41126
_____	Kit, SurePath™ 1 kit = 25 vials/cervex-brush combi #	45031
_____	Kit, ThinPrep® 1 kit = 25 vials/brooms #	12587
_____	Kit, ThinPrep® 1 kit = 25 vials/brushes/spatulas #	40624
_____	Kit, UroVysion FISH #	41440
_____	Mailers, Two Slide PK/24	16474

MEDIA/CULTURE KITS

_____	Anaerobic Transport Media PK/10	16426
_____	Kit, Aptima Combo 2 - Urine Collection PK/50	28908
_____	Kit, Aptima Combo 2 - Unisex Swab PK/50	28907
_____	Kit, BD AFFIRM VPIII PK/10 #	29624
_____	Kit, Aptima, Combo 2 (Thin Prep)	42711
_____	Kit, Chlamydia Collection DFA	16578
_____	Kit, QuantiFERON-TB (High Altitude)Gold In-Tube	46261
_____	Kit, QuantiFERON-TB (Standard) Gold In-Tube	45112
_____	Kit, Stool Transport (Unpreserved) PK/40	40910
_____	Kit, Timed Stool Collection, if needed also order (#28077, Container, Fecal Fat: 1 for 48hr or 2 for 72hr)	44192
_____	Media, Enteric Transport C&S (Cary-Blair) PK/20	29799
_____	Media, Regan-Lowe Bordetella Pertussis	24962
_____	Media, UTM, Swab (M4) PK/50	12884
_____	RPMI Media	24967
_____	Swab, Liquid Aimes (eSwab) PK/50	45877
_____	Swab, Nasopharyngeal Ultrafine PK/100	41334
_____	UBIT Breath Test	51124

Issuing unit of measure is "each" unless specified in the item description as PK or CS.

#—ITEM IS GROUND SHIPPED ONLY

QUANTITY

SHIPPING CONTAINERS

_____	Box, Shipping (10" x 10" x 8")	16546
_____	Infectious Shipping Container, Medium	19814

COLLECTION CONTAINERS: **NOT FOR TRANSPORT**

****WILL BE DISCARDED IF SUBMITTED****

_____	Urine, 24 Hour (Graduated, w/ Lid Handle)	16419
-------	---	--------------

SPECIMEN CONTAINERS

_____	Container, Fecal Fat, Stool, 1L	28077
_____	Pyruvate (Perchloric Acid) #	16567
_____	Rack, Transport Tube (10-Hole) EACH	50779
_____	Rack, Transport Tube (10 Hole) CS/400	10772
_____	COVER, for 10-Hole Rack	28094
_____	Rack, Transport Tube (20 Hole)	42769
_____	COVER, for 20-Hole Rack	42770
_____	Serotonin (Ascorbic Acid)	16568
_____	Tube, Amber Standard Transport PK/100	13654
_____	Tube, Microtainer, Pediatric Lead K2EDTA	24347
_____	Tube, Standard Transport, CS/1000	15824
_____	Tube, Transport, Trace Metals Free PK/100	43116
_____	Tube, 4 mL Sterile Transport PK/100	43115
_____	Tube, 15 mL Conical (Sterile) PK/100	11282
_____	Tube, 50 mL Conical (Sterile) PK/50	29582
_____	Vacutainer™, 8.5 mL Yellow (ACD Sol A)	16380
_____	Vacutainer™, 6 mL Royal Blue (EDTA)	16313
_____	Vacutainer™, 6 mL Royal Blue (Plain) (Metal Free)	16359
_____	Vacutainer™, 8.3 mL SPS for Microbiology	24964
_____	VIP/Glucagon/PTHRP	49662

MISC/ADDITIONAL SUPPLIES

_____	Boric Acid (1 gram tablets) Urine Preserve PK/100	53481
-------	---	--------------

ARUP provides nonstandard containers and supplies ONLY for the collection and transport of lab specimens intended for ARUP Laboratories.

